



Printing Solutions

Omnicell Printing Solutions Streamline Your Labeling Process

Rising medication costs. Growing concerns around patient safety.
Increased pressures on your staff.

As the healthcare landscape continues to evolve, it is more important than ever to support safe, accurate medication delivery to patients while realizing efficiencies where possible.

Omnicell's printing solutions are designed to meet all printing and labeling solutions of your institutional pharmacy, from unit-dose foil printing to license plate labels and label application as well as color label printing, including medication images.



Fast, Flexible and Easy Medication Labeling

License plate labeling system

The license plate labeling system perfectly meets the needs of the long-term care pharmacy that requires barcode-scannable labels. Whether you need horizontal labels for standard punch cards, vertical labels for unit dose boxes special labels for syringes, refrigerated items, compounds and more—this system gives you the flexibility to create it all.

- Design options for identification of narcotics, day of the week, delivery times and more
- Peel-and-present or batch printing options
- Printable daily production report



Reduce medication waste with more packaging options

Unit dose foil labeling system

Foil lidding containing all pertinent drug information and barcodes, helping you improve efficiency and control costs while meeting the needs of your customers.

- Accommodate diverse packaging needs including unit-dose box, short-cycle dispensing, reclaimable blister cards, eKits and automated dispensing cabinets
- Reduce waste by pre-packing expensive medications into reclaimable blister cards (as allowed by state regulations)
- Track log numbers and expiration dates using barcode technology, with option of printable daily production report



Assurance at a glance

Create high-resolution, full-color labels for fast, accurate visual verification with the PillVue™ color labeling system.

- Link images to manufacturer's barcode data, reducing potential for data entry errors
- Create images for new medications at any time
- Color-code labels for quick visual identification of expiration dates or time pass
- Differentiate yourself by making your customers' lives easier



Labels can be color-coded for fast identification of expiration dates.

The printing solutions backbone: EasyLabel Software

Omniceil's reliable and easy-to-use Windows-based software supports all styles of Omnicell packaging used in long-term care pharmacies. The software's unique multi-view printing (MVP) feature enables simultaneous management and printing of several label types to streamline the packaging process.

- Unit-dose Foil
- License Plate Labels
- PillVue Color Labels
- Label Conveyor Applicator



EasyLabel software enables compliance and efficiencies through the following capabilities:

- FDA 21 CFR Part II compliant
- 2D barcode-ready
- Windows XP or Windows 7 (recommended) compatible
- 16-digit, RSS/GS1 barcode and RFID capable
- Three licensing options available:
 - Terminal: supports multiple printers from multiple workstations among multiple locations
 - Network: supports multiple printers from multiple workstations within a single location
 - Single user: supports three printers per workstation (when connected via USB, serial or parallel port, the maximum distance between the computer and any printer is 15 feet)

Technical requirements for Omnicell's printing solutions

	EasyLabel software	License plate system	Unit dose	PillVue color labeling
Windows 7 or higher	✓	✓	✓	✓
CD ROM drive	✓			
1 USB port for each printer, the scanner, and the pill-imaging camera (if needed)	✓	✓	✓	✓
Network drop (if communicating via Ethernet)	✓			
Standard printer for report printing	✓	✓	✓	✓

Break through bottlenecks in pre-pack with the Label Conveyor Applicator (LCA) System

The LCA System automates the labeling process for your pre-pack operation while delivering reduced labor time and increased pre-pack efficiency. Using a high-speed label applicator system, the label conveyor applicator with the LCA scan functionality is used immediately after the packing process to automate the application of the pre-pack label. It prints, positions and sticks the label, guaranteeing perfect placement every time.

Key benefits include:

- Ability to export tech production report to track individual technician productivity
- Ability to export production report to track inventory



Production report

NDC NUMBER	QTY	DRUG NAME	STRENGTH	LOT NUMBER	TECH_INT	SCANNED	PRODUCTS PER MIN
65862059805	14	TAMSULOSIN	0.4MG CAP	TMSA15447A	MA	98	1,075 ppm
68382023016	14	FENOFIBRATE	145MG TAB	1053074	MA	58	2,041 ppm
68382023016	14	LANSOPRAZOLE	30MG CAP	406086	MA	45	3,640 ppm
51991077290	28	MAGNESIUM OXIDE	400MG TAB	48154	MA	34	1,152 ppm
57896053101	28	MULTIVITAMIN W/ MINERALS	TAB	531F01	MA	32	12,315 ppm
00904198860	60	MAPAP ES	500MG tab	47523	MA	30	0,727 ppm

Tech report

Tech: MA	QTY	DRUG NAME	STRENGTH	PRODUCTS PER MIN
14	FENOFIBRATE	145MG TAB	58	
14	LANOPRAZOLE DR	30MG TAB	98	
14	TAMSULOSIN	0.4MG CAP	98	
28	MAGNESIUM OXIDE	400MG TAB	34	
28	MULTIVITAMIN W/ MINERALS	TAB	32	
60	MAPAP ES	500MG TAB	30	
	TOTAL		297	

Visit Omniceil.com/Medication-Adherence to learn more today.

For Your Pharmacy. For Your Patients.

Your partner for optimized outcomes.

Through industry-leading medication packaging and automation solutions, Omnicell empowers pharmacies to work more efficiently, to improve medication safety and satisfaction for their employees, customers and patients.

