


Omnicell

Integrated Automation Solutions




Omnicell

Founded in 1992, Omnicell is a leading provider of technologically advanced automation that enables healthcare facilities to acquire, manage, dispense, and deliver medications and supplies more effectively.

Omnicell is exclusively focused on medication and supply dispensing automation. That means we understand the workflow and problems faced by pharmacy, nursing, materials management, and IT. With our newly designed fourth-generation Omnicell G4 platform, we are making it easier than ever before for our customers to reduce costs, keep pace with ever-shifting regulatory requirements and safeguard the patient.

Our passion for what we do is evident by the unsurpassed customer loyalty we have earned. Omnicell offers first-rate customer support, a comprehensive training and education program, and a range of professional consulting services.

Whether you are implementing medication and supply automation for the first time or planning to upgrade or expand current systems, count on Omnicell to be your trusted partner.



Omnicell: Innovating for the Future

For nearly two decades, Omnicell has been laser-focused on one thing: developing automation to improve efficiency, safety, and regulatory compliance in hospitals.

Our new G4 automated medication and supply dispensing platform builds on this foundation by offering new features and functionality that help hospitals raise the bar in each of these areas. Products across the G4 platform operate with a single, unified database. The dispensing and reconciliation of medication from controlled substance storage to all patient care areas of the hospital, including bedside, operating rooms, and procedural areas helps provide a closed-loop system.

In addition to this seamless integration, the G4 platform introduces a larger, clearer display and redesigned interface that improve user experience with the automated dispensing cabinet.

Recent Omnicell Innovations

Omnicell's fourth-generation G4 platform brings remarkable innovation to 11 different medication and supply automation products and solutions, including:

- A new **G4 Console** with a larger, high-resolution display, 2-D bar code capability for accurate medication verification, a state-of-the-art biometric fingerprint ID system, an integrated medication label printer, an **8-bin Locking Drawer**, and the latest **Omnicell 15.0** software.
- **Savvy™ Mobile Medication System** provides safe, secure transportation of medications from the automated dispensing cabinets to the bedside.
- New **Controlled Substance Management** system provides perpetual inventory management, compliance with regulatory standards, and increased efficiency for the pharmacy.
- Redesigned **Anesthesia Workstation™** keeps anesthesia supplies and medications securely stored yet easy to access during surgery.
- Updated **OptiFlex™ 11.0** supply system to decrease inventory costs while reducing stock-outs by easing the process of replenishment.
- New **PandoraVIA® Cabinet Management Module** automatically manages stockouts and PAR inventory levels as well as helps identify potential diversion.
- New partnerships with three leading companies in healthcare: **CardinalASSIST®** expands options in cabinet restock, **RxScan®** for bar code check and verification on all doses leaving the pharmacy, and **Sentry Data Systems** to reduce medication costs to hospitals that participate in the 340B program.




www.klasresearch.com

Best in KLAS Award: OmniRx: 2010

KLAS Product Leadership Award:
OmniRx: 2006–2010
OmniLinkRx: 2009–2010
WorkflowRx and Carousel: 2010



Central Pharmacy



The battle cry to *do more with less* seems never-ending, but Omnicell can help your pharmacy make this a reality by alleviating the burden of compliance documentation and giving your staff time for higher priority patient safety issues.

Our **Carousel**, **Bar Code Packager**, and **WorkflowRx™** software deliver a powerful, automated, inventory management system to increase productivity and reduce costs. Without adding staff, you can achieve greater control and tracking of medications, improve dispensing accuracy, and enable dynamic reordering of all pharmacy stock, including 340B medications.

OmniLinkRx™ manages and prioritizes the vast number of physician orders that are written every day, expediting delivery from all patient care areas to the pharmacy — and eliminating interruptions caused by lost or misplaced paperwork.

Omicell's completely new **Controlled Substance Management (CSM) system** improves management of narcotics. CSM uses the same database as your Omnicell dispensing cabinets. With CSM, Omnicell helps provide closed-loop tracking of medication during dispensing and reconciliation to ensure a complete audit trail throughout the medication distribution process.

SafetyStock® bar code verification software with **RxScan** provides barcode verification of restock orders in the Central Pharmacy.

Omicell has partnered with Cardinal Health to offer the **CardinalASSIST** program — an automated process to replenish drugs in specific units of measure for stocking in Omnicell cabinets.

Our partnership with **Sentry Data Systems** helps hospitals reduce medication costs.

"From receiving to dispensing, the WorkflowRx solution is a capability that allows us to increase our control over medication inventory levels. Pharmacy operations were streamlined by 25% to 50% and efficiency based on user-defined IOR scheduled rules increased by 50% to 75%."

—David Grant, RPh, MBA,
Director of Pharmacy
Chambersburg Hospital



Perioperative

Anesthesia Workstation

Interoperability

OptiFlex SS

Tissue Tracking

OmniRx Automated Medication and Supply System



"The benefits from materials running smoother in the OR and improved efficiencies far outweigh the investment we made in both time and cost... I would ask to have [Omnicell] installed in any hospital I worked for."

—Geoff Hibbert,
OR Manager
Patewood Memorial Hospital

In the perioperative areas, Omnicell offers flexible solutions for managing medications, tissue, and supplies. By eliminating manual processes, these systems dramatically increase charge capture while preventing delays caused by out-of-stock items.

The **Anesthesia Workstation**—now with a more compact design—is an automated dispensing cabinet for the operating room that helps you meet medication security regulations while meeting the needs of anesthesia providers.

The **Omnicell Tissue Center** helps hospitals comply with Joint Commission and AORN standards to improve patient safety by enabling efficient, electronic documentation and management of tissue specimens.

The **OptiFlex™ SS (Surgical Services) system** integrates physician preference cards with perpetual inventory management to streamline supply inventory management—improving efficiency for clinical and supply chain staff alike.

OptiFlex CL

Interoperability

OmniRx
Automated
Medication
and Supply
System

SinglePointe

Hanger CARES

Anywhere RN

Omniceil helps you get the right medication or supply to the patient on time and in the most efficient manner. We offer innovations you won't find anywhere else to streamline workflow and prevent medication errors.

Omniceil's modular, flexible **Automated Dispensing Cabinets** keep your automation assets on the cutting edge. The G4 cabinet console offers the new **Touch & Go™** biometric system for increased efficiency and security, an integrated medication label printer, a larger high-resolution display, and Microsoft® **Windows 7** operating system and the latest Intel® Core™ multi-core processor which improves performance and assures long-term support.

The **OmniDispenser™** single-dose dispensing system eliminates tedious narcotics count-backs.

SinglePointe™ is our unique software that enables virtually all medications for each patient to be stocked in the automated dispensing cabinet.

Anywhere RN™ gives nurses the freedom to remotely access the Omnicell cabinet and manage medications from anywhere in a patient care area – reducing trips to the cabinet and minimizing interruptions that could lead to medication errors.

OptiFlex™ MS allows you to manage medical-surgical supplies on open shelves, closed cabinets, or a combination of the two. **OptiFlex CL™** is designed for the unique needs of specialty areas such as the cath lab or endoscopy. Both systems greatly improve inventory management and charge capture.

The unique **CARES® vendor-managed inventory program**, offered through Omnicell's partnership with Hanger Prosthetics & Orthotics is a secure, automated, cost-effective solution for dispensing orthopedic soft goods and select durable medical equipment.



"Implementation of new automated patient-specific management software has improved patient safety by decreasing the medication error rate associated with the patient transfer process, improved operational efficiencies for pharmacy, and provided for more consistent nursing practice."

—Sherri Ramsey,
Director of Pharmacy,
University of Tennessee
Medical Center


Mobile

Anywhere RN

Savvy
Mobile
Medication
System

"Anywhere RN unchains nurses from the cabinet and gives them more time to review medication orders prior to accessing medications. This flexible and adaptable platform improves productivity by giving our nurses control over their own workflow."

—Joseph Alessandrini, R.Ph.,
Director of Pharmacy,
South Jersey Healthcare




Omniceil's mobile systems add another level of versatility and safety to your medication distribution process by creating a versatile link between the automated dispensing cabinet and the patient.

The Omnicell **Savvy™ Mobile Medication System** allows nurses to safely transport medications from the automated dispensing cabinet to the patient rooms. The system reduces the number of trips back and forth to the cabinet, and it improves patient safety by securing medications and enabling verification at the bedside.

Savvy has our **Anywhere RN™** remote medication management software built in, so nurses can remain in the patient room to queue up and waste medications remotely.

Business Analytics Interoperability



Omnicell offers analytics tools to help hospitals gain valuable insight from the data generated by their dispensing systems.

PandoraVIA® is our latest business analytics tool that delivers on-demand, easy-to-follow reports analyzing your medication use and supply chain processes. PandoraVIA helps you prioritize action items that keep you in compliance, lower costs, and improve patient care.


In addition to PandoraVIA, Omnicell offers other web-based analytics products and customized reporting as well.

Interoperable solutions integrate products, provide drug accountability, and increase patient safety and efficiency.

Omnicell's interoperable, standards-based systems effectively complement other vendor solutions in the hospital inpatient setting. Our systems are able to share clinical information between disparate systems to eliminate redundant data entry and enhance reporting capabilities.

"The data analytics from Pandora greatly enhance the reporting tools on our Omnicell cabinets to help identify drug diversion and more efficiently meet compliance regulations. With Pandora, we are able to identify diverters much quicker, hopefully before it becomes an issue of patient safety."

—Patty Womack
Director of Pharmacy



Our Obsession is Your Solution.

1201 Charleston Road
Mountain View, CA 94043
Phone: 800.850.6664
Fax: 650.251.6266
www.omnicell.com

3661 Bur Wood Drive
Waukegan, IL 60085
Phone: 847.596.3400
Fax: 847.596.3402

443 Donelson Pike
Suite 200
Nashville, TN 37214
Phone: 615.564.7299
Fax: 615.564.7297

10 Victor Square
Suite 250
Scotts Valley, CA 95066
Phone: 866.429.8900
Fax: 831.458.2123

