

Pharmacy Drawer Options

Omnicell offers a wide variety of drawer options for storing, controlling, and accurately tracking medications and pharmacy supplies. Pharmacy drawers integrate into the Omnicell® automated dispensing cabinets and Anesthesia Workstation™ to create customized systems for specific clinical needs.

Flexible Cabinet Configurations

Omnicell drawers accommodate many types of medications, including oral solids and liquids, syringes, ampules, vials, IV solutions, and PCA medications. The drawers can be configured within the cabinet in a variety of ways.

- Drawers are available in varying degrees of security and varying bin sizes.
- Pharmacy drawers and medical supply areas can be combined in the same cabinet.

Cabinet Type	Drawer Capacity
OmniTT™ tabletop cabinet	Up to 3 drawers
OmniRX® cabinet	Up to 13 drawers*
Each cell of the one-cell, two-cell, and three-cell cabinets	Up to 18 drawers*

*An OmniDispenser module takes up 6 drawers.

OmniDispenser™ Single-Dose Dispensing System

The OmniDispenser stores medications in a separate, secured location in the cabinet and dispenses each dose individually, so nurses only have access to the specific quantities requested. The system automatically tracks every dose, virtually eliminating time-consuming countbacks to confirm remaining inventories of controlled substances. This is a significant time saver for nurses that gives them more time for patient care.

Key Features

- Virtually eliminates controlled substance countbacks and related discrepancies.
- Push-button technology helps ensure accurate restocks.
- Facilitates improved compliance with Joint Commission and other regulatory standards.

Guiding Lights

Omnicell's unique guiding lights technology visually directs clinicians to the drawer and then to the appropriate bin for the medication they had selected. This feature reduces the possibility of nurses taking the wrong item.





FlexBin™ Single-Dose Drawer

With this drawer, access to a medication is restricted to the item and quantity chosen by the clinician, providing a high level of security. No other doses are available. There is no need to confirm the remaining quantity or conduct end-of-shift counts, which results in fewer discrepancies.

Key Features

- Guiding lights illuminate drawer location and individual bins
- Single-deep drawer maximizes inventory in a minimum amount of space
- Three bin sizes available to accommodate a variety of medication types
- Bins can be used as single-dose or multi-dose
- Flexible drawer configurations: 96 small bins, 48 medium bins, and 48 small/large combinations

Locking Lid Drawer

Locking bins, suitable for controlled substances, provide a high level of security by restricting access to one preselected medication at a time. An audio alarm sounds if a clinician attempts to access a non-preselected compartment and notification is recorded and sent to pharmacy. The tamper-proof metal lids of each bin further ensure security.

Key Features

- Guiding lights illuminate drawer location and individual bins
- Single-deep and double-deep drawers
- Four bin sizes: 6 bins, 8 bins, 12 bins, and 24 bins
- Tamper-proof, locked metal lids

Lighted Matrix Drawer

The Lighted Matrix Drawer provides a lower level of security but maximum flexibility. Lights frame the bin and guide the user to the preselected medication, reducing the possibility of errors.

Key Features

- Guiding lights illuminate drawer location and individual bins
- Flexible configuration with adjustable bin dividers
- Capacity up to 24 unique line items

Matrix Drawer

The Matrix Drawer is a high-capacity, low-security drawer that is suitable for large quantities of non-controlled medications.

Key Features

- Adjustable bin dividers with capacity up to 96 unique items

Omnicell Medication/Supply Drawer

The Omnicell Medication/Supply Drawer provides a low level of security in width-configurable open bins. Placing this drawer within a pharmacy zone created in the cabinet can further control access to pharmacy items.

Key Features

- Guiding lights illuminate drawer location and individual bins
- Easy-to-use green reorder buttons
- Adjustable bin dividers
- Capacity up to 24 unique line items



Omnicell, the Omnicell design mark, and OmniRx are registered trademarks. Anesthesia Workstation, FlexBin, OmniDispenser, and OmniTT are trademarks of Omnicell, Inc. in the United States and internationally. All other trademarks and trade names are the property of their respective owners. ©2015 All rights reserved. PN709 Rev.B 06/15