

Omniceil Anesthesia-Rx

PROVIDERS CHOICE FOR A FAST, INTUITIVE WORKFLOW IN THE OR

- Complements existing OR workflow
- Supports fast, reliable access to medication
- Improves synergy between pharmacy and OR
- Helps increase safety and compliance



Designed specifically for the OR

Omniceil's Anesthesia-Rx™ automated medication cart complements existing operating room workflow to provide quick access to medication with just a few touches of the screen. Anesthesia-Rx, designed in collaboration with anesthesia providers, enhances communication between pharmacy and the OR while helping to support safety and regulatory compliance.



Complements workflow with fast, reliable access

Anesthesia-Rx* helps expedite medication delivery to the patient by offering a variety of features to help users quickly find and dispense medications. An intuitive, easy-to-use interface and touch screen reduces time spent interacting with the cart, allowing providers to focus on the patient. Secure, large-capacity drawers are designed so medications, including controlled substances, are easily accessible, while high-capacity drawer lids accommodate larger items such as Propofol vials.



Anesthesia-Rx eliminates the need for manual documentation through an electronic data collection process that improves charge capture, waste capture and controlled substance reconciliation. In the pharmacy, the workload for checking, charging and restocking is minimized, making documentation easier to manage and significantly reducing the amount of paperwork.

Improves synergy between OR and pharmacy

As medications are stocked, dispensed, returned and wasted, pharmacy is able to view reports to efficiently track and optimize inventory levels to help reduce stock outs. Expiration Date Tracking provides the ability to remove expired medications from the cart to avoid medication errors.

Anesthesia-Rx improves security through drug utilization accountability. The cart captures transactions, making controlled substance utilization and waste verification more efficient. The ability to monitor drug inventory and document individual access helps reduce controlled substance discrepancies and prevent and detect diversion, while supporting regulatory compliance.

Supports safety and compliance

Anesthesia-Rx features the only integrated color label printer on the market, which reduces the time anesthesia providers spend preparing medications, and helps to improve regulatory compliance and patient safety. When a medication is dispensed, a color label can be automatically printed that meets regulatory guidelines. Color coding allows for easy identification to reduce medication identity confusion.

Lockout prevention options allow providers to quickly access medications in the event of an emergency. These include a Manual Lock that will not lock once a case has begun, a Manual Override Button designed for emergency access, a traditional key lock for back up access, and an uninterrupted power supply to maintain access during a power outage.

To learn more about Anesthesia-Rx, please contact your Omnicell representative or visit www.Omnicell.com.

* Anesthesia-Rx is a Class 1 Medical Device and adheres to strict quality standards as defined by the FDA.

Omnicell and the Omnicell design mark are registered trademarks of Omnicell, Inc. in the United States and internationally. Anesthesia-Rx is a trademark of Aesynt Incorporated. All other trademarks and trade names are the property of their respective owners.
©2016 All rights reserved. MM-ASRX-1 / Revision: 1 / DCR-632 / Effective Date: 17 June 2016