

Savvy Mobile Medication Workstation



Mobile. Reliable. Secure.
A smarter way to manage medications.

Savvy Mobile Medication Workstation

Improve Nurse Efficiency to Increase Patient Satisfaction

Patients are more satisfied with their care when their nurses spend more time with them. The Savvy™ mobile medication workstation creates a safer, simpler workflow that supports less rushed interactions between your nurses and patients.

Are your nurses spending enough time with patients?

30%

Nurses spend nearly 30% of their time engaged in activities not directly tied to patient care.¹⁻²

Are interruptions causing medication errors?

13%

Every time a nurse is interrupted, the chance of a clinical error occurring increases nearly 13%.³

Is nurse fatigue contributing to medication errors?

5

Nurses often walk up to 5 miles per 12-hour shift.⁴



Savvy Solution

Savvy wirelessly links in real time to the automated dispensing cabinets (ADCs) and includes Anywhere RN™ remote medication management software, allowing nurses to remotely manage medications in quieter areas where they are less likely to be interrupted.

With Savvy, nurses can make a single trip to the ADC for multiple patients, reducing steps and leaving more time for patient care.

When interruptions do occur during a med pass, nurses have peace of mind knowing that medications are secured in patient-assigned locking drawers.

Streamlined Medication Management

- Facilitates faster medication order processing
- Shortens waiting times at the ADC
- Allows remote requests for ordering, returning, and wasting of medication
- Ensures 24/7 availability with hot-swappable batteries — eliminates plug-in recharging

Heightened Safety and Security

- Addresses Institute for Safe Medication Practices (ISMP) recommendations for safe transport of medications
- Increases medication accountability between the ADC and patient bedsides
- Eliminates manual drawer labeling—patient drawers are automatically assigned
- Increases patient safety with independently locking drawers—simply scan the patient ID bar code to unlock patient-specific drawer
- Guiding lights technology guides nurses to correct patient drawer

33%

Savvy Advantage

Nurses achieved a 33% decrease in ADC medication removal time.⁵

77%

Savvy Advantage

77% of nurses were able to double-check medications due to time savings with Anywhere RN.⁵



Seamless Integration Assures Efficiency and Accuracy

- Savvy wirelessly communicates in real-time with the Omnicell ADC, which is integrated with the hospital information systems
- Shared database with Omnicell medication dispensing system streamlines data entry, reducing the potential for errors

Savvy Nurses, Satisfied Patients

The Savvy mobile medication workstation increases medication safety while cutting down on steps and interruptions that disrupt patient care. With reimbursements becoming increasingly tied to patient satisfaction,⁶ shouldn't you get Savvy about medication management?

¹ Lewis N. Process management tools for nursing lacking. *Information Week*, April 1, 2010. ² Hendrich A, Chow M, Skierczynski B, Lu Z. A 36-hospital time and motion study: how do medical-surgical nurses spend their time? *Permanente Journal*. Summer 2008;12:25-34. ³ Westbrook J, Woods A, Rob M, Dunsmuir W, Day R. Association of interruptions with an increased risk and severity of medication administration errors. *Archives of Internal Medicine*. April 26, 2010. 170(8):683-690. ⁴ How far do nurses walk? *MedSurg Nursing*, August 2006. ⁵ Data on file. ⁶ Adamy J. U.S. ties hospital payments to making patients happy. *Wall Street Journal*. October 14, 2012. Accessed online 2/5/13.

Savvy Mobile Medication Workstation

Specifications

Display

Input	12 VDC
Resolution	1280 x 1024 (maximum)
Type	TFT active matrix with touch screen; privacy filter
Size	17" diagonal

Integrated Computer

Processor	Intel SP9300 CPU; 2.26 GHz C2D, 6M L2 cache, ultra low power
Operating System	Windows 7
Memory	4 GB, 800 MHz DDR2 SO-DIMM
Data Storage	80 GB SATA SSD
Video	Integrated SXGA 1280 x 1024 32-bit color
Ports	2 external USB, 1 internal USB, 1 Ethernet, 1 DVI
Wireless	Intel Centrino Advanced – N6205, 802.11 a/b/g/n; WiFi wireless adapter

Power System

Input Voltage	90 – 264 VAC @ 3.5 A maximum at 115 VAC, 1.6 A maximum at 230 VAC
Output Voltage	120 VAC @ 60 watts maximum
Internal Charger	4 amps maximum
Battery Type	Lithium-ion, hot swappable, up to 18 hours runtime (2 batteries)
Short-Circuit Protection	AC 4 A fuse; 12 VDC 7 A fuse; motor 10 A fuse; inverter 5 A resettable fuse
Status Indicator	On-screen display; LED indicators at base; LED indicators on each battery case

Savvy External Battery Charger

Dimensions	16.75" W x 11" H x 7.75" D
Number of Charging Bays	2
Recharge Time	5-6 hours from fully depleted condition
Input Voltage	100-240 VAC @ 3.0 A maximum at 115 VAC, 1.5 A maximum at 230 VAC
Approvals	TUV, IEC60601-1, FCC

General

Base Dimensions	16" x 16"
Work Surface Area	17" x 12"
Work Surface Height	36" to 46"
Device Height	Adjustable height between 55"-65" (including monitor)
Construction	Cast aluminum and plastic injection
Wheels	4 swivel (2 locking) casters, 5" diameter, polyurethane tread
Approvals	CE/CSA, UL60601-1, FCC
Drawers	Up to 12 locking drawers with Omnicell guiding lights technology; 3" and 6" deep drawer options

To learn more about the Savvy mobile medication workstation, please contact your Omnicell representative or visit www.Omnicell.com.

See a brief video demonstration of the Savvy mobile medication workstation.



Omnicell and the Omnicell design mark are registered trademarks. Anywhere RN and Savvy are trademarks of Omnicell, Inc. in the United States and internationally. All other trademarks and trade names are the property of their respective owners. ©1999-2013 All rights reserved. PN714 Rev-D