

## Omniceil XT Automated Dispensing Cabinets and Drawers

No other medication dispensing system offers more cabinet choices to meet the needs of acute and post-acute care sites. The scalable Omnicell® XT Automated Dispensing System can easily be expanded with additional cabinets, and drawers can be added or changed on-site as clinical needs evolve. Medications and supplies can be combined in a single cabinet, promoting convenience and saving space.



# Omnicell XT Automated Dispensing Cabinets



## Three-Cell Cabinet

### Dimensions:

- Height: 77.5"
- Width: 76.5"
- Depth: 27.0"



## Two-Cell Cabinet

### Dimensions:

- Height: 77.5"
- Width: 51.5"
- Depth: 27.0"



## One-Cell Cabinet

### Dimensions:

- Height: 77.5"
- Width: 26.5"
- Depth: 27.0"



## Half-Height Cabinet

### Dimensions:

- Height: 53.3"
- Width: 26.5"
- Depth: 27.0"



## Quarter-Height Cabinet

### Dimensions:


- Height: 25.4"
- Width: 26.5"
- Depth: 27.0"



# Omniceil XT Cabinet Drawers



10-Bin Drawer



18-Bin Drawer



36-Bin Drawer



27-Bin Drawer



6-Bin Double-Deep Drawer



24-Bin Double-Deep Drawer



48-Bin Drawer

## Metal Locking Lid Drawers

Durable metal locking lids keep medications secure. All drawers include Guiding Lights.

**10-Bin:** Larger bin fits up to three 60 mL PCA units.

**Bin Type: 5 Extended-Wide**

- Height: 1.7"
- Width: 6.6"
- Length: 3.4"

**Bin Type: 5 Long**

- Height: 1.7"
- Width: 13.8"
- Length: 3.4"

**18-Bin:** Securely stores a variety of items, including prefilled syringes.

**Bin Type: Extended**

- Height: 1.7"
- Width: 6.6"
- Length: 2.6"

**36-Bin:** Securely stores oral solids, vials, and small par level items.

**Bin Type: Standard**

- Height: 1.7"
- Width: 2.9"
- Length: 2.6"

**27-Bin:** Larger bin accommodates prefilled syringes and other larger items.

**Bin Type: 18 Standard**

- Height: 1.7"
- Width: 2.9"
- Length: 2.6"

**Bin Type: 9 Extended**

- Height: 1.7"
- Width: 6.6"
- Length: 2.6"

**6-Bin:** Largest locked bin in the industry—fits IV bags and other bulky items.

**Bin Type: Double-Deep**

- Height: 4.0"
- Width: 10.1"
- Length: 6.3"

## Open Configurable Drawers

Now with greater configurability and brighter Guiding Lights.

**24-Bin:** Double-deep drawer can be configured to include up to 24 bins.

**48-Bin:** Drawer can be configured to include up to 48 bins.

# Omnicell XT Auxiliary Cabinets



## Three-Cell Auxiliary Cabinet

### Dimensions:

- Height: 77.5"
- Width: 76.5"
- Depth: 27.0"

## Two-Cell Auxiliary Cabinet

### Dimensions:

- Height: 77.5"
- Width: 51.5"
- Depth: 27.0"

## One-Cell Auxiliary Cabinet

### Dimensions:

- Height: 77.5"
- Width: 26.5"
- Depth: 27.0"

Are you ready to make medication management easier?  
Learn more about the XT Series at [www.omnicell.com/xt](http://www.omnicell.com/xt).