

# How can you achieve convenient access to medications and compliance in the OR?

## Omnicell XT Series Anesthesia Workstation

Securing medications in the operating room is about more than patient safety. It is about protecting the entire surgical team, the pharmacy, and the reputation of the hospital. We understand the high-pressure, life and death situations that anesthesiologists experience every day. We also understand the regulatory requirements the pharmacy must account for.

Omnicell's Anesthesia Workstation™ gives anesthesia providers instant access to medications while providing pharmacy with tighter control in the operating room.

### Access Anesthesia Providers Want


Fast, easy access to medications for anesthesia providers:

- **Instant Access** – Unique breakaway bar enables emergency access to medications
- **Fast and Secure** – Advanced biometric ID provides fast access with military grade security
- **Convenient** – High capacity anesthesia cabinet helps prevent stockouts that slow case turnaround
- **Efficient** – Automates reconciliation of controlled substances and eliminates countbacks

### Security Pharmacy Needs

Pharmacy gains tighter control of medication usage in the OR:

- **Compliant** – Secures controlled substances at all times and provides greater inventory oversight
- **Efficient** – Tracks inventory usage and expiration dates
- **Secure** – Employs Windows 10 Enterprise Data Protection
- **Safe** – Enables accurate restocking with SafetyStock® bar code confirmation option



## Maximize Patient Safety While Maintaining Compliance

Anesthesia Workstation securely stores medications and supplies for a full day of cases. It automatically tracks inventory used, cutting hours of manual documentation.

For added patient safety, Anesthesia Workstation supports the ability to generate color-coded medication labels, which helps reduce the three most common medication errors made in the OR: vial/ampoule swaps, mislabeling, and syringe swaps. Medication label printers support compliance with the following regulatory standards and guidelines:

- Joint Commission NPSG.03.04.01
- ISO 26825 Regulations
- ASTM 4774 Guidelines
- ASA Standards

## Bring It All Together With Unity

As part of the Unity enterprise platform, the Anesthesia Workstation integrates with Omnicell automation systems deployed in the central pharmacy, nursing units, and other locations. This unified approach helps to reduce IT complexity, minimize system maintenance, and prevent errors. It also enables more comprehensive inventory visibility and diversion monitoring insights.



Learn more about how the XT Anesthesia Workstation can help you meet the demands of the surgical suite by visiting [www.omnicell.com/aws](http://www.omnicell.com/aws).